



STACK

A.E.D Checklist

FORM ID:	517607
STATUS:	Completed
CREATED ON:	01/17/2025 2:38 AM (UTC-05:00)
COMPLETED ON:	01/17/2025 2:39 AM (UTC-05:00)
ELAPSED TIME:	0 Hour(s) 0 Minute(s)
COMPLETED BY:	David Bedillion
GPS LOCATION:	Not Specified
ASSET:	Philips DEN-AED
FACILITY:	Denver HQ Office 1700 Broadway Ste. 1750 Denver, CO 80202
LOCATION:	Other - Operations
LOCATION DETAILS:	Outside RM 117A
NOTES:	

General Information

Pass

Is the AED readily visible and not obstructed from view?

Pass

Is the AED unit located in the designated area with proper signage?

Pass

Is the AED unit clean and undamaged?

Pass

If the AED in an alarmed cabinet, is the alarm functioning?

Pass

Check AED Supplies

Pass

Is the AED supplied with a CPR supply kit (i.e. disposable razor, scissors, pocket mask, gloves, etc.)?

Pass

Does the kit contain a disposable razor (to remove chest hair that interferes with the application of pads)?

Pass

Does the kit contain scissors (to cut difficult to remove garments)?

Pass

Does the kit contain a CPR pocket mask (note expiration date)?

Pass

Does the kit contain a package of gloves?

Pass

Check AED Supplies - Expiration Dates (CPR Pocket Mask)

QUESTIONS WITH FORMATTED NUMBERS & TEXT:

Expiration Date - CPR Pocket Mask

-

AED Status Indicator

Pass

Is the green light visible on the unit? If the light is red, notify the AED program coordinator. If the indicator light is absent check the battery.

Pass

Defibrillator Electrode Pads

Pass

Are there light blue and white wires running from the sealed defibrillator pad package?

Pass

Are there two sets of pads available?

Pass

Do the packages appear to be sealed, cleaned, and undamaged?

Pass

Find the expiration date on the pads, are they within their expiration date?

Pass

Defibrillator Electrode Pads - Expiration Dates

QUESTIONS WITH FORMATTED NUMBERS & TEXT:

Electrode Pads Set 1 - Expiration Date

-

Electrode Pads Set 2 - Expiration Date

-

Batteries

Pass

When opening the unit, does the display indicate a green blinking light in the top right corner above the green power button?

Pass

Does the unit have a spare battery (not required)? If the unit has a spare battery, follow steps for testing using the AED inspection guide.

Pass

Batteries - Expiration Dates

QUESTIONS WITH FORMATTED NUMBERS & TEXT:

Batteries - Expiration Dates

-